




UNI 31

FEATURES

- Low service costs
- Long life
- Short term storage of blood and plasma
- Digital controller with potential free contact (which makes it possible to connect it to an external alarm)



UNI

PERFECT FOR RELIABLE
STORING IN LOW TEMPERATURE



UNI 11



UNI 21



UNI 41



UNI 51

	UNI 11	UNI 21	UNI 31	UNI 41	UNI 51
Amb. Temp. range deg. (°F)	+50 to +86	+50 to +86	+50 to +86	+50 to +86	+50 to +86
Gross volume, cu. ft	4.6	8.1	10.7	12.8	14.9
Net volume, cu.ft	4.6	8.1	10.7	12.8	14.9
Temperature range inside unit (°F)	+14 to -50	+14 to -50	+14 to -50	+14 to -50	+14 to -50
Refrigerant: R290	Yes	Yes	Yes	Yes	Yes
Input power (watt):	308	310	330	370	385
Energy consumption / 24 h (kwh)	1.5	1.9	1.9	2.5	2.6
DIMENSIONS					
External dimensions (wxdxh)	28.5x25.8x34.1	41.5x25.8x34.1	51.4x25.8x34.1	59.3x25.8x34.1	67.1x25.8x34.1
Internal dimensions (wxdxh)	20.5x17.7x25.6	33.5x17.7x25.6	35.4x17.7x25.6	43.3x17.7x25.6	51.2x17.7x25.6
Dimension A	28.5	41.5	51.4	59.3	67.1
Dimension B	25.8	25.8	25.8	25.8	25.8
Dimension C	34.1	34.1	34.1	34.1	34.1
Dimension D	30.7	30.7	30.7	30.7	30.7
Dimension E	20.5	33.5	43.3	51.2	59.1
Dimension F	12.6	25.6	35.4	43.3	51.2
Dimension G	17.7	17.7	17.7	17.7	17.7
Dimension H	25.6	25.6	25.6	25.6	25.6
Dimension J	15.7	15.7	15.7	15.7	15.7
Dimension L	33.5	33.5	33.5	33.5	33.5
Insulation, mm	3.9	3.9	3.9	3.9	3.9
ACCESSORIES					
Dividers	Optional	Optional	Optional	Optional	Optional
Max. number of baskets (option)	1	2	3	4	4
Castors	Optional	Optional	Optional	Optional	Optional
Castors w. brakes	Optional	Optional	Optional	Optional	Optional
Lock	Included	Included	Included	Included	Included
Digital controller w. acoustic alarm	Included	Included	Included	Included	Included
OTHER DATA					
Outer cabinet	White	White	White	White	White
Inner cabinet	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Type of cooling	Static	Static	Static	Static	Static
Type of defrost	Manual	Manual	Manual	Manual	Manual
Power 110volt / 60 hz	Yes	Yes	Yes	Yes	Yes
Potential free contact	Included	Included	Included	Included	Included
Thermostat	Digital	Digital	Digital	Digital	Digital
Net weight, lbs	101	130	147	163	176
Gross weight, lbs	108	136	156	172	187
SHIPPING					
20' cont. Capacity	42	30	24	20	18
40' cont. Capacity	90	64	50	46	40
Standard Trailer capacity	102	74	60	52	44